Epi Navigator Cath Flex

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INDICATIONS:

The Epi Navigator Cath Flex epidural catheter is used in percutaneous epidural lysis and/or adhesiolysis procedures of the spinal canal, which is reached through the dedicated access needle.

This procedure can be associated with the possibility of treating roots and ganglia exiting the liberated spinal foramen with Pulsed Radiofrequency (PRF) - (prior to spinal cord stimulation - SCS)

The Epi Navigator epidural catheter is also equipped with a thermocouple for the real-time monitoring of the temperature reached by the target nervous tissue and an infusion lumen for localised pharmacological treatment (prior to administration of analgesic drugs).

This surgical procedure is indicated for patients with symptoms of:

- FBSS
- Spinal Stenosis
- Chronic Low Back Pain
- Radiculopathy
- Postherpetic Neuralgia

Caution: Patients with pathologies other than those listed are treatable, with responsibility and effectiveness attributable to the operator.

Caution: For correct use of the device, carefully read the following information.

CHARACTERISTICS:

Epi Navigator Cath Flex is distally equipped with a **5-mm** active electrode which enables the supply of the current required for Pulsed Radiofrequency (PRF) treatment and an internal lumen for the infusion of drugs through the proximal hemostatic valve (infusion is also possible with stylet inserted).

The ergonomic handle allows the operator to direct the distal tip and better position the active electrode.

Once the desired foramen has been reached, it is possible to perform stimulation tests on the sensory and motor nerve fibres.

Caution: Positioning of the tip must always be carried out under scopic control.

During PRF treatment, a temperature sensor continuously monitors the thermal condition of the active zone.

Following mechanical adhesiolysis and Pulsed Radiofrequency treatment, it is possible to leave the catheter inside the patient (by connecting it to the extension provided in the package) to proceed with repeated infusions of analgesics in the post-operative phase.

This device therefore makes it possible to perform:

- Radiofrequency treatments for neuromodulation of nervous tissue;
- Single or repeated injections into the spinal cavities where pharmacological treatment is required (NB: the device allows the infusion of drugs for 48 hours after surgery in the operating theatre thanks to the dedicated extension provided in the package);
- Tests for bone marrow stimulation before implantation of a spinal cord stimulation system;
- Tests for application of morphine/cortisone based solutions before implantation of a pumping system;
- Stimulation of peripheral nerves to detect the origin of pain;
- Percutaneous removal of adhesions (percutaneous epidurolysis and/or adhesiolysis of the spinal canal).

SIDE EFFECTS:

- Infections (apply the best possible aseptic conditions and avoid any sources of infection. It is advised to perform the usual perioperative antibiotic prophylaxis);
- Hemorrhages and hematomas (minimise risks by operating with special precautions and analysing the patient's coagulogram);
- Thermal damage to structures following an incorrect setting of the instrument;
- Allergic reactions due to medical substances injected (indications for the application of used medications must be observed, taking into account the associated side effects);
- Mechanical/structural lesions (minimisable with accurate preparation and careful execution of the operation);
- Functional complications due to dislocation, breaking or shearing of the catheter;
- Do not use the device on women who are pregnant or exposed to pregnancy;

Caution: Single-use, non-resusable device.

WARNINGS:

- **Check that compatible generators are set in "PULSED RADIOFREQUENCY (PRF)" mode;**
- **NEVER set instruments in "THERMAL LESION" or "CONTINUOUS RADIOFREQUENCY" mode;**
- Check that the data set on the generator are as follows (to be considered as standard parameters for use):
- **Temperature: 42 °C Frequency: 2Hz Duration of stimulation impulse: 20 ms Total duration of treatment: 240 sec;**
- Tissue impedance is obtained automatically from the generator based on the tissue adjacent to the electrode;
- Observe the Generator instructions for use and warnings.

INTRODUCTION OF THE SYSTEM:

1. Attach the dispersion plate to the patient (for example, along the calf...);
2. Using fluoroscopic images, place the dedicated needle at the introduction point of the catheter into the epidural space;
3. Extract the stylet leaving the sleeve in place and insert the catheter through this;
4. Carefully advance the catheter into the epidural space, performing adhesiolysis manoeuvres if the spinal tract presents stenotic regions;

Caution: Keep the impedance parameter monitored. A value less than 200Ω could imply that the electrode has been placed in the subdural region.

If you wish to curve the distal portion of the catheter, extract the stylet up to the black mark and move the lever on the handle according to the direction to be taken. Advance the catheter into the desired area by pushing it through the needle and checking its movement through radiographic images;

Caution: Do not reinsert the stylet into the catheter when it is angled.

To reinsert the stylet it is necessary to straighten the catheter using the handle, remove the catheter from the patient, lay it on a sterile surface and advance the stylet through the dedicated lumen, taking care not to pierce the catheter.

6. When the desired target has been reached, connect the electrocatheter to its stimulation device, using the cable supplied with the product.

Connection cables: **Ref: VE-100-XX** where **XX** identifies the first two letters of the Generator to which they are destined (BA: Baylis – Kimberly Clark; CO: Cosman; NE: Neurotherm; RA: Radionics; ST: Stryker);

Caution: Observe the stimulation device manufacturer's indications.

7. Following the operation, and after the various lysis and PRF treatments on the ganglia concerned, it is possible to inject pharmacological solutions from the tip of the catheter through the dedicated two-way valve;

8. If, instead, you wish to leave the catheter in place for infusions in the **48 hours following** the operation, proceed as follows:

Completely extract the stylet from the body of the catheter; cut the catheter along the indicated red mark; remove the introducer needle taking care not to move the catheter; connect the irrigation extension to the end of the catheter using a one-way purple hemostatic valve; proceed with drug infusion using a Luer Lock syringe; fix the external portion of the catheter with Steri-Strips and precise dressing.

Caution: To avoid infection during the post-operative phase, apply sterile dressings and replace frequently.

REMOVAL OF THE SYSTEM:

At the end of the procedure, extract the catheter and cover the introduction point with an antiseptic cream and sterile dressing.

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