
Epi Navigator Needle

IT	AGO INTRODUTTORE PER EPI NAVIGATOR CATH Istruzioni per l'uso
EN	INTRODUCER NEEDLE FOR EPI NAVIGATOR CATH Instructions for use
DE	EINFÜHRKANÜLE FÜR EPI NAVIGATOR CATH Bedienungsanweisung
ES	AGUJA INTRODUCTORA PARA EPI NAVIGATOR CATH Instrucciones de uso

Rev. 00/2024

AMS Group S.r.l. s.u.
Via Europa, 12
35020 San Pietro Viminario (PD) - Italia
Tel. +39 0429 719378 - Fax +39 0429 760135

CE 0051

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35020 San Pietro Viminario (PD) - Italia

Tel. +39 0429 719378 - Fax +39 0429 760135

EPI NAVIGATOR NEEDLE**Introducer needle for Epi Navigator Cath and Epi Navigator Cath Flex****INDICATION:**

Introducer Needle for insertion into the epidural space of Epi Navigator Cath and Epi Navigator Cath Flex, the electrocatheter used for epidural diagnostic and acute and chronic pain therapy.

CHARACTERISTIC:

Device with epidural cannula tip, provided with luer lock and inner stylet.

TECHNICAL DATA:

- Diameter: 14G
- Length: 100mm
- Pipe: AISI 304 e 302
- Disposable device. Not reuse.
- Ethylene Oxide (ETO) sterilization.
- Double packaging in medical paper bag and laminated plastic film.
- Sterility is guaranteed only to undamaged packages.
- Do not use if opened or damaged packaging.
- Device to be used only by medical staff with specialist knowledge and only in health care.
- Do not disperse in the environment.
- Store away from direct heat and protected from sunlight.
- Read carefully and follow Epi Navigator Cath and Epi Navigator Cath Flex instructions for use and warnings.

COLLATERAL EFFECTS:

- Infections (apply the best sterile conditions in order to avoid possible sources of infection. We recommend carrying out common pre-operation antibiotic prophylaxis).
- Hemorrhage and hematoma (minimize risks by using particular caution and execute coagulogram tests);
- Allergic reactions due to injected medical substances (instructions for use of medicines used must be respected, taking into consideration relative collateral effects).
- Structural mechanical lesions (these can be minimized with careful preparation for the operation's execution – the majority of introduction methods are standardized).

Y: Do not executed this treatment on women who are pregnant or risk pregnancy.

Y: Disposable material: not to be re-sterilized. Risk of infections.

INTRODUCTION AND REMOVAL SYSTEM:

- Verify packaging and device integrity before use;
- the best sterile conditions in order to avoid possible sources of infection;
- Remove the protective cap from the needle insert the introducer into the epidura space;
- Remove the inner stylet only when properly positioned into the epidural space and before proceeding with the introduction of the electrocatheter;
- Pay attention to any movement of the catheter in the opposite direction of the needle tip, to avoid possible brackage of the catheter on the cutting edge of the needle;
- Insert the catheter and proceed carefully;
- Do not retract the catheter once inserted into the needle, in order to avoid accidental cuts of the same catheter;
- After treatment, remove the needle and dispose in appropriate containers;
- To avoid infection during catheter use, apply sterile medication.
- Read carefully and follow Epi Navigator Cath and Epi Navigator Cath Flex instructions for use and warnings.

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