GRAPHIC SYMBOLS USED IN LABEL			
Symbol	Description	Symbol	Description
	Do not re-use – Single use		Use by date
	Batch code		Sterilized using ethylene oxide
	Catalogue number		Consult instructions for use
	Do not use if package is damaged		Manufacturer
	Do not re-sterilize		EC certification according MDD 93/42
	Non-pyrogenic		Date of manufacture

AMS Group S.r.l. denies all responsibility for any damage caused to the patient and/or medical staff if the system is used anomalously, improperly or due to mistaken movements and unqualified staff.

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4. Insert the video-guided catheter into the introducer previously placed during the access procedure.
5. Inject a small quantity of sterile saline solution.
6. Push the video-guided catheter forward under fluoroscopic guidance and apply a light wash if necessary, then proceed to observe the area on the video monitor.

Warning: "UP", "DOWN", "LEFT" and "RIGHT" marks near the orientation levers correspond to the real movements of the catheter, however they may not fully correspond to the endoscopic images on video screen (especially if the images are not aligned horizontally as explained at paragraph 2).

Warning: the "OUT" liquid collection line may be connected to the relevant extension line and/or to the liquid collection bag.

The extension line is provided with **Resascope**, while the collection fluid bag is contained in the SK-ST-10 procedural kit.

Warning: the lines marked with "TOOL1" and "TOOL2" are intended for the insertion of IESS® surgical instruments, such as the RESAFLEX Quantic Molecular Resonance¹ terminal and the RESALOON balloon catheter, and other mechanical micro-instruments. Follow the relevant instructions for use of each IESS® accessory.

Warning: in the soft tip version of the **Resascope**, before inserting the required surgical instruments, completely remove the stiffening stylet previously inserted in the "TOOL2" channel (stylet use is limited to the access into the epidural space phase).

7. If necessary, use the line marked with "IN" to administer medications or other medical solutions.

Warning: do not administer medications before extracting the optical fiber flexible endoscope or any other powered instrument and closing the relevant valves of the dedicated channel.

8. Remove the **Resascope** video-guided catheter and the introducer sheath.

Warning: before inserting the device into the patient, make sure that the system works properly: puller wires, valves, "LOCK" and lens with its focus.

Warning: Verify bending limits on the instruction manual of the optical fiber endoscope. Excessive bending may cause permanent damage to the endoscope.

The orientation lock systems of the tip are activated by sliding forward the specific sliders marked with "LOCK". To return to the standard position slide back the "LOCK" inserted in the opposite direction.

Warning: Do not operate the moving levers when the "LOCK" system is activated.

Warning: Do not use if package is damaged. Sterile and non-pyrogenic. Single use.

Warning: The tip of the video-guided catheter **Resascope** may be configured into 2 different versions (soft tip and rigid tip), specifically in the version with rigid tip operate the moving levers with care as the pulling wires are subject to higher traction.

INSTRUCTIONS FOR USE

Warning: for the correct positioning of the sheath, required to insert the video-guided catheter in the epidural space, follow the instructions for use provided with the access kit (percutaneous introduction method into the sacral hiatus or e posterolateral hiatus).

1. Connect a 10cc or 20cc prefilled saline syringe to the stopcock identified with the "IN" mark (the saline will help to enlarge the epidural space and ensure a better endoscopic view of the epidural space).

Warning: the package contains 2 extension lines to increase, if necessary, the length of the infusion and reflux ways; a collection bag for drained fluids may also be connected to the extension lines.

Warning: saline must be infused slowly into the epidural space to avoid blood pressure peaks.

2. Insert the sterile optical fiber flexible endoscope into the video-guided catheter through the line marked with "OPTIC" and slide the endoscope forward up to the catheter tip. Verify that the distal tip of the flexible endoscope is not obstructed. When the correct position is reached, block the optical fiber through at the entrance "OPTIC" channel and verify that the images on the screen are properly positioned (this is usually done by framing some written signs).

Warning: Do not apply excessive strength when closing the valve as this could damage the flexible endoscope.

3. Using prefilled syringe, fill the video-guided catheter with saline and remove any air bubble.

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INSTRUCTIONS FOR USE

RESASCOPE EPIDURAL STEERABLE VIDEO-GUIDED CATHETER

INDICATIONS

The video-guided catheter **Resascope**, used with a flexible optical fiber endoscope, allows to:

- Explore the epidural space and evaluate its anatomical and pathological conditions
- Infuse washing solutions and medications into the epidural space
- Operate diagnostic and surgical tools to treat and investigate several epidural pathologies

This surgical procedure is indicated for patients with symptoms related to:

- FBSS
- Spinal Stenosis
- Chronic Low Back Pain
- Radiculopathy

Warning: Patients with other pathologies, other than those listed, can be treated under the operator's responsibility and effectiveness.

CONTRAINDICATIONS

The **Resascope** video-guided catheter is contraindicated when, at the physician's discretion, its use can threaten the patient's health, due to known allergies to its components.

Warning: spinal anaesthesia is contraindicated for procedures involving minimally invasive surgical approach of the epidural space (spinal endoscopy).

Do not use distilled water solutions!

Possible complications: as in every surgical procedure possible complications may occur during postoperative phase.

Specifications of this procedure may include: headache, retinal or intracranial hemorrhage, dural puncture, meningoenitis, discitis.

Possible side effects: possible side effects may include pain at the insertion site and paresthesia/dysesthesia in the perioperative phase.

GENERAL INSTRUCTIONS

The tip of the **Resascope** video-guided catheter can be flexed by using the two levers on the handle.

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AMS Group S.r.l. s.u.